

Akicare Technology Limited

Email: kevin@akicare.com.cn WhatsApp:+8613622363037

Reusable Biopsy Guides List

Brand	Support Probe Model	Probe Type	Model	Materials	Angle	Support Needle Size
ALPINION	EN3-10	Endocavity	JEM-027	Stainless Steel	Single Angle	16-18G
ALPINION	L3-12, L3-12T, L3-12i, L3-12H	Linear	JSM-143	Stainless Steel	Three Angles	11-23G
ALPINION	SC1-4H, C1-6, C1-6T, C1-6i, C1-6CT, SC1-6, SC1-6H	Convex	JSM-159	Stainless Steel	Three Angles	11-23G
ALPINION	L3-12X, L8-17X	Linear	JSM-199	Stainless Steel	Three Angles	11-23G
ALPINION	E3-10, E3-10H	Endocavity	JEM-030	Stainless Steel	Single Angle	16-18G
ALPINION	EC3-10, EV3-10, EC3-10H, EC3-10X, EV3-10H, EC3-10T, EV3-10T, EV3-10X	Endocavity	JEM-033	Stainless Steel	Single Angle	16-18G
ALPINION	EV2-11H, EC2-11H	Endocavity	JEM-095	Stainless Steel	Single Angle	16-18G
ALPINION	VE3-10H	Endocavity	JEM-094	Stainless Steel	Single Angle	16-18G
BK medical	8570, 8670	Linear	JSM-286	Stainless Steel	Three Angles	14-22G
BK medical	8819, E10C4, 9019	Endocavity	JEM-030	Stainless Steel	Single Angle	16-18G
BK medical	8823, 8802, 8803, 8820e, 8830, 6C2, 9040, 6C2s, 9023, 9C2, 9002	Microconvex	JSM-005	Stainless Steel	Single Angle	14-22G
BK medical	8848, 8808, E14CL4b, 9048, 8658S & 8658T	Transvaginal	JSM-283	Stainless Steel	Single Angle	14-18G
BK medical	8811, 8822, 10L2w, 9022, 13L4w, 9011	Linear	JSM-233	Stainless Steel	Three Angles	11-23G
BK medical	8808e, 8818, E14C4t, 9018, E11C3b, 9008	Transvaginal	JSM-377	Stainless Steel	Single Angle	14-18G
BK medical	8808e, 8818, E14C4t, 9018, E11C3b, 9008	Rectal	JEM-058	Stainless Steel	Single Angle	16-18G
BK medical	8667, E14C4	Rectal	JEM-077	Stainless Steel	Single Angle	16-18G
Canon	PVF-641VT, PVM-651VT, PVN-661VT, PVQ-641V, PVT-661VT, PVG- 630v	Endocavity	JEM-018	Stainless Steel	Single Angle	16-18G
Canon	PST-25AT, PST-25BT, PST-25ST, PST-30BT, PST-30SBT, PSU-30BT	Phased	JSM-085	Stainless Steel	Two Angles	11-23G
Canon	PVT-375AT, PVT-375AX, PVT-375BT, PVT-375ST, PVU-375BT, PVT-375SC	Convex	JSM-009	Stainless Steel	Two Angles	11-23G
Canon	PLT-604AT, PLT-704AT, PLT-704SBT, PLT-704ST, PLU-704BT	Linear	JSM-017	Stainless Steel	Two Angles	11-23G
Canon	PVF-375AT, PVK-357AT, PVM-375AT, PVQ-375A, PVN-375AT	Convex	JSM-086	Stainless Steel	Three Angles	11-23G

Canon	PVT-382BT, PVU-382BT	Microconvex	JSM-027	Stainless Steel	Two Angles	11-23G
Canon	PLT-805AT, PLT-1005BT, PLU-1005BT	Linear	JSM-068	Stainless Steel	Two Angles	11-23G
Canon	PLK-703AT, PLM-703AT, PLF-703AT, PLN-703AT,	Linear	JSM-095	Stainless Steel	Two Angles	11-23G
Canon	PSK-25AT, PSM-25AT, PSM-30BT, PSM-37AT, PSM-37CT, PSM-50AT, PSM-70AT, PSN-25AT, PSN-37AT, PSN-37CT, PSN-50AT, PSN-70AT	Phased	JSM-098	Stainless Steel	Two Angles	14-22G
Canon	PLT-1204AT, PLT-1204BT, PLU-1204BT	Linear	JSM-108	Stainless Steel	Single Angle	11-23G
Canon	PVF-381MT, PVM-381AT, PVM-662AT, PVM-672AT	Microconvex	JSM-169	Stainless Steel	Two Angles	14-22G
Canon	PVT-781VT, PVU-781VT, PVT-781VTE, PVU-781VTE	Endocavity	JEM-019	Stainless Steel	Single Angle	16-18G
Canon	PLG-805S	Linear	JSM-202	Stainless Steel	Two Angles	11-23G
Canon	PVT-681MVL, PVU-681MVL	Endocavity	JEM-005	Stainless Steel	Single Angle	16-18G
Canon	PVT-681MV	Endocavity	JEM-032	Stainless Steel	Single Angle	16-18G
Canon	PVT-770RT, PVM-740RT	Endocavity	JEM-061	Stainless Steel	Single Angle	16-18G
Canon	PVT-350BTP, PVU-350BTP, 6CP1	Convex	JSM-261	Stainless Steel	Four Angles	18-19G
Canon	PVT-674BT, PVU-674BT, PVQ-674	Convex	JSM-227	Stainless Steel	Two Angles	11-23G
Canon	PVT-712BT, PVU-712BT	Microconvex	JSM-267	Stainless Steel	Three Angles	11-23G
Canon	PVL-715RT, PVL-715RS, PVL-715RST	Transvaginal	JSM-269	Stainless Steel	Single Angle	18G
Canon	PLT-705BT	Linear	JSM-155	Stainless Steel	Three Angles	11-23G
Canon	PVI-475BX (i8CX1)	Convex	JSM-209	Stainless Steel	Five Angles	11-23G
Canon	PLI-2004BX (i24LX8)	Linear	JSM-218	Stainless Steel	Five Angles	11-23G
Canon	PLI-1205BX (i18LX5)	Linear	JSM-256	Stainless Steel	Five Angles	11-23G
Canon	PSI-30BX	Phased	JSM-304	Stainless Steel	Free Angle	11-23G
Canon	PVI-475BT (i8C1), PVT-475BT	Convex	JSM-314	Stainless Steel	Five Angles	11-23G
Canon	PVT-482BT(8MC1), PVI-482BX(i8MCX1)	Microconvex	JSM-335	Stainless Steel	Five Angles	16G/18G
Canon	PVU-621VST(10C4), PVU-781VST(11C4)	Endocavity	JEM-069	Stainless Steel	Single Angle	16-18G
Canon	PVM-621VT, PVG-601V, PVF-621VT, PVF-651VT	Endocavity	JEM-085	Stainless Steel	Single Angle	16-18G
Canon	PVU-770ST	Endocavity	JEM-089	Stainless Steel	Single Angle	16-18G
Canon	PLI-705BX(i11LX3)	Linear	JSM-379	Stainless Steel	Five Angles	11-23G
Canon	PVT-482BT(8MC1), PVI-482BX(i8MCX1)	Microconvex	JSM-384	Stainless Steel	Five Angles	11-23G
CHISON	D3C60L, C3-A, C3-E, C3-V	Convex	JSM-241	Stainless Steel	Three Angles	11-23G
CHISON	D3C20L, MC3-A, MC3-E, MC3-V	Microconvex	JSM-325	Stainless Steel	Three Angles	11-23G
CHISON	32C60LA	Convex	JSM-300	Stainless Steel	Three Angles	11-23G

CHISON	D7L40L, L7-E, L7-V, D12L40L, L12-E, L12-V	Linear	JSM-196	Stainless Steel	Three Angles	11-23G
CHISON	D6C12L, D7C10L, V6-A, V6- E, V7-E, E6-V, E7W-V	Endocavity	JEM-069	Stainless Steel	Single Angle	16-18G
EDAN	P5-1XQ	Phased	JSM-026	Stainless Steel	Two Angles	14-23G
EDAN	C361-1, C361-2, C362, C363-1	Convex	JSM-042	Stainless Steel	Three Angles	14-23G
EDAN	C5-1Q, C5-2Q	Convex	JSM-288	Stainless Steel	Three Angles	14-23G
EDAN	E611, E611-1, E611-2, E613, E612UB, E8-4D, E8-4Q	Endocavity	JEM-015	Stainless Steel	Single Angle	16-18G
EDAN	MC8-4Q, MC8-4D	Microconvex	JSM-263	Stainless Steel	Three Angles	14-23G
EDAN	C5-2Q, C5-2D, C5-2XQ	Convex	JSM-234	Stainless Steel	Three Angles	14-23G
EDAN	L12-5Q(212623011, 212634011), L12-	Linear	JSM-240	Stainless Steel	Three Angles	14-23G
EDAN	L17-7HQ, L17-7HD, L10-4Q, L10-4D, L12-5Q, L12-5D	Linear	JSM-273	Stainless Steel	Four Angles	14-23G
EDAN	P5-1Q	Phased	JSM-276	Stainless Steel	Two Angles	14-23G
EDAN	MC9-3TQ, MC9-3TD	Microconvex	JSM-184	Stainless Steel	Two Angles	14-23G
EDAN	E10-3BQ	Endocavity	JEM-066	Stainless Steel	Single Angle	16-18G
EDAN	E10-3HQ	Endocavity	JEM-067	Stainless Steel	Single Angle	16-18G
EDAN	L12-5HQ	Linear	JSM-301	Stainless Steel	Three Angles	14-23G
EDAN	MC8-4Q, MC8-4D	Microconvex	JSM-263	Stainless Steel	Three Angles	14-23G
EDAN	L12-5WQ	Linear	JSM-343	Stainless Steel	Three Angles	14-23G
EDAN	ECL12-3Q	Transvaginal	JSM-347	Stainless Steel	Single Angle	18G
EDAN	L12-5HQ	Linear	JSM-412	Stainless Steel	Three Angles	14-23G
EDAN	BC4-160H	Convex	JSM-437	Stainless Steel	Three Angles	11-23G
EDAN	E10-3MQ	Endocavity	JEM-091	Stainless Steel	Single Angle	16-18G
Esaote	EC123, EC1123, SE3123, SE3133, E3-12	Endocavity	JEM-020	Stainless Steel	Single Angle	16-18G
Esaote	EC1123, SE3123, SE3133, E3- 12	Endocavity	JEM-012	Stainless Steel	Single Angle	16-18G
Esaote	CA421, CA430, CA430E, CA431	Convex	JSM-002	Stainless Steel	Two Angles	14-22G
Esaote	LA522, LA523, LA532	Linear	JSM-012	Stainless Steel	Three Angles	11-23G
Esaote	LA332, LA332E	Linear	JSM-107	Stainless Steel	Three Angles	11-23G
Esaote	E8-5	Endocavity	JEM-028	Stainless Steel	Single Angle	16-18G
Esaote	CA541	Convex	JSM-053	Stainless Steel	Three Angles	11-23G
Esaote	CA621, CA631	Convex	JSM-106	Stainless Steel	Two Angles	14-22G
Esaote	CA123	Microconvex	JSM-029	Stainless Steel	Two Angles	11-23G
Esaote	TRT33, TLC3-13	Transvaginal	JSM-193	Stainless Steel	Single Angle	18G

Esaote	SL1543, L4-15, LA533, AL2443	Linear	JSM-197	Stainless Steel	Two Angles	11-23G
Esaote	AC2541, C1-8	Convex	JSM-198	Stainless Steel	Three Angles	11-23G
Esaote	AL2442, L3-11	Linear	JSM-253	Stainless Steel	Three Angles	11-23G
Esaote	LA435, LA424, SL2325	Linear	JSM-113	Stainless Steel	Three Angles	11-23G
Esaote	C5-2R13, SC3121	Microconvex	JSM-124	Stainless Steel	Free Angle	18-20G
Esaote	SC3123, CA1123	Microconvex	JSM-125	Stainless Steel	Two Angles	11-23G
Esaote	mC3-11	Microconvex	JSM-125	Stainless Steel	Two Angles	11-23G
Esaote	PA250, P1-5	Phased	JSM-340	Stainless Steel	Free Angle	14-22G
Esaote	CX1-8	Convex	JSM-421	Stainless Steel	Three Angles	11-23G
Esaote	LMX3-15	Linear	JSM-422	Stainless Steel	Three Angles	11-23G
Esaote	LMX4-20	Linear	JSM-423	Stainless Steel	Three Angles	11-23G
FUJIFILM Healthcare	EUP-U531	Rectal	JEM-016	Stainless Steel	Single Angle	16-18G
FUJIFILM Healthcare	EUP-V33W, EUP-V53W, C41V	Endocavity	JEM-034	Stainless Steel	Single Angle	16-18G
FUJIFILM Healthcare	EUP-CC531, EUP-CC531s, CC41R, CC41R1	Endocavity	JEM-034	Stainless Steel	Single Angle	16-18G
FUJIFILM Healthcare	EUP-C314T, EUP-C314G, EUP-C324T, EUP-C514, EUP-	Convex	JSM-079	Stainless Steel	Four Angles	11-23G
FUJIFILM Healthcare	EUP-C514	Convex	JSM-079	Stainless Steel	Single Angle	11-23G
FUJIFILM Healthcare	EUP-L54MA	Linear	JSM-371	Stainless Steel	Three Angles	11-23G
FUJIFILM Healthcare	EUP-L74M, L55	Linear	JSM-018	Stainless Steel	Three Angles	11-23G
FUJIFILM Healthcare	UST-981-5, UST-984-5	Endocavity	JEM-068	Stainless Steel	Single Angle	16-18G
FUJIFILM Healthcare	EUP-B512	Microconvex	JSM-030	Stainless Steel	Two Angles	11-23G
FUJIFILM Healthcare	UST-979-3.5, UST-989(D)-3.5, UST-990-5, UST-9114-3.5, UST-9115-5, UST-9119, UST-9123, UST-9126, UST-9127, UST-9130	Convex	JSM-033	Stainless Steel	Single Angle	11-23G
FUJIFILM Healthcare	UST-5412	Linear	JSM-052	Stainless Steel	Three Angles	11-23G
FUJIFILM Healthcare	UST-5410, UST-5410H, UST-5413	Linear	JSM-052	Stainless Steel	Three Angles	11-23G
FUJIFILM Healthcare	UST-5284-2.5, UST-5285-3.5, UST-5286-2.5, UST-5287-3.5, UST-5292-5, UST-5294-5, UST-5297, UST-5299, UST-	Phased	JSM-102	Stainless Steel	Three Angles	11-23G
FUJIFILM Healthcare	EUP-C715, C25	Convex	JSM-037	Stainless Steel	Four Angles	11-23G
FUJIFILM Healthcare	UST-9118, UST-9124	Endocavity	JEM-059	Stainless Steel	Single Angle	16-18G
FUJIFILM Healthcare	EUP-C532, C42	Microconvex	JSM-114	Stainless Steel	Three Angles	11-23G
FUJIFILM Healthcare	UST-9133, C22K	Microconvex	JSM-055	Stainless Steel	Two Angles	11-23G
FUJIFILM Healthcare	UST-52105	Phased	JSM-133	Stainless Steel	Three Angles	11-23G

FUJIFILM Healthcare	EUP-U533, C41L47RP	Transvaginal	JSM-138	Stainless Steel	Single Angle	18G
FUJIFILM Healthcare	UST-944B-3.5, UST-9102U-3.5	Microconvex	JSM-178	Stainless Steel	Single Angle	11-23G
FUJIFILM Healthcare	EUP-C311	Microconvex	JSM-187	Stainless Steel	Two Angles	14-22G
FUJIFILM Healthcare	UST-9128, UST-978-3.5, UST-9103-5, UST-9121-3.5	Microconvex	JSM-200	Stainless Steel	Single Angle	11-23G
FUJIFILM Healthcare	EUP-L52, L34, UST-5419	Linear	JSM-201	Stainless Steel	Three Angles	11-23G
FUJIFILM Healthcare	EUP-B712, C22P	Microconvex	JSM-211	Stainless Steel	Two Angles	16G/18G
FUJIFILM Healthcare	EUP-V73W, C41V1	Endocavity	JEM-050	Stainless Steel	Single Angle	16-18G
FUJIFILM Healthcare	C251, C35, C252, C253	Convex	JSM-070	Stainless Steel	Four Angles	11-23G
FUJIFILM Healthcare	UST-5415, L441	Linear	JSM-151	Stainless Steel	Three Angles	11-23G
FUJIFILM Healthcare	EUP-S70, S21	Phased	JSM-242	Stainless Steel	Free Angle	18-20G
FUJIFILM Healthcare	UST-9101-7.5	Convex	JSM-213	Stainless Steel	Single Angle	11-23G
FUJIFILM Healthcare	EUP-L73s, L44	Linear	JSM-215	Stainless Steel	Single Angle	11-23G
FUJIFILM Healthcare	EUP-B514	Convex	JSM-014	Stainless Steel	Three Angles	16-18G
FUJIFILM Healthcare	C41B	Endocavity	JEM-056	Stainless Steel	Single Angle	16-18G
FUJIFILM Healthcare	EUP-B715, C25P	Convex	JSM-121	Stainless Steel	Three Angles	18-19G
FUJIFILM Healthcare	EUP-L75, L64, UST-5417	Linear	JSM-054	Stainless Steel	Three Angles	11-23G
FUJIFILM Healthcare	S211, S12	Phased	JSM-212	Stainless Steel	Free Angle	11-23G
FUJIFILM Healthcare	UST-9147	Convex	JSM-067	Stainless Steel	Four Angles	14-22G
FUJIFILM Healthcare	EUP-L65	Linear	JSM-140	Stainless Steel	Three Angles	11-23G
FUJIFILM Healthcare	UST-670P-5, UST-675P, UST-676P, C41RP	Rectal	JEM-051	Stainless Steel	Single Angle	16-18G
FUJIFILM Healthcare	UST-9120, UST-987-7.5, UST-9104-5, C42K	Microconvex	JSM-266	Stainless Steel	Single Angle	11-23G
FUJIFILM Healthcare	UST-672-5/7.5, UST-678	Transvaginal	JSM-134	Stainless Steel	Single Angle	14/16G, 18/20G
FUJIFILM Healthcare	UST-5268P-5, UST-52114P, S31KP	Phased	JSM-128	Stainless Steel	Single Angle	11-23G
FUJIFILM Healthcare	L35	Microconvex	JSM-335	Stainless Steel	Five Angles	16G/18G
FUJIFILM Healthcare	CL4416R, CL4416R1	Transvaginal	JSM-395	Stainless Steel	Single Angle	18G
FUJIFILM Sonosite	L25, L25e, L25x, L25xp	Linear	JSM-076	Stainless Steel	Free Angle	18-20G
FUJIFILM Sonosite	HFL38xp, L38xi, L38xp, HFL38xi, (L12-3)	Linear	JSM-117	Stainless Steel	Three Angles	11-23G
FUJIFILM Sonosite	C60, C60e, C60n, C60x	Convex	JSM-156	Stainless Steel	Single Angle	11-23G
FUJIFILM Sonosite	C8x	Rectal	JEM-001	Stainless Steel	Single Angle	16-18G
FUJIFILM Sonosite	ICT, ICTe, ICTn, ICTx, ICTxp	Endocavity	JEM-024	Stainless Steel	Single Angle	16-18G
FUJIFILM Sonosite	HFL38, HFL38x, L38, L38e, L38x	Linear	JSM-165	Stainless Steel	Single Angle	11-23G

FUJIFILM Sonosite	P21x, P21xp	Phased	JSM-220	Stainless Steel	Two Angles	11-23G
FUJIFILM Sonosite	HFL50x, HFL50xp	Linear	JSM-280	Stainless Steel	Single Angle	11-23G
FUJIFILM Sonosite	C60xp	Convex	JSM-349	Stainless Steel	Three Angles	11-23G
FUJIFILM Sonosite	C60xf, C60xi, rC60xi, (C5-1)	Convex	JSM-399	Stainless Steel	Three Angles	11-23G
GE Healthcare	E8C, E8C-RS, E8CS, E8CS-RS, E8CS-RS, E8CS-SC, IC5-9H, IC9B-RS	Endocavity	JEM-037	Stainless Steel	Single Angle	16-18G
GE Healthcare	E8C(Voluson), IC5-9-D	Endocavity	JEM-057	Stainless Steel	Single Angle	16-18G
GE Healthcare	4C, 4C-A, 4C-D, 4C-RS, 4C-SC	Convex	JSM-003	Stainless Steel	Three Angles	11-23G
GE Healthcare	11L, 12L, 11L-D, 12L-RS, 12L-SC, L6-12-RS	Linear	JSM-010	Stainless Steel	Three Angles	11-23G
GE Healthcare	L4-12t-RS, L12n-RS, L12n-SC	Linear	JSM-010	Stainless Steel	Three Angles	11-23G
GE Healthcare	8L, 8L-RS	Linear	JSM-077	Stainless Steel	Three Angles	11-23G
GE Healthcare	RIC5-9-D, RIC6-12-D, RIC5-9W-RS, RIC5-9A-RS, RIC5-9-RS, RIC5-9H, RIC5-9	Endocavity	JEM-063	Stainless Steel	Single Angle	16-18G
GE Healthcare	3.5C (348C), 548C, 3.5CS, AC2-5	Convex	JSM-097	Stainless Steel	Single Angle	11-23G
GE Healthcare	C1-5-D, C1-5-RS	Convex	JSM-025	Stainless Steel	Three Angles	11-23G
GE Healthcare	S1-5,S1-5-D	Phased	JSM-100	Stainless Steel	Two Angles	11-23G
GE Healthcare	12L(LA39)	Linear	JSM-101	Stainless Steel	Single Angle	11-23G
GE Healthcare	ML6-15, ML6-15-RS, ML6-15-D	Linear	JSM-036	Stainless Steel	Three Angles	11-23G
GE Healthcare	IC5-9	Endocavity	JEM-011	Stainless Steel	Single Angle	16-18G
GE Healthcare	9L, 9L-D, 9L-RS	Linear	JSM-048	Stainless Steel	Three Angles	11-23G
GE Healthcare	M5S, M5S-D	Phased	JSM-035	Stainless Steel	Three Angles	11-23G
GE Healthcare	3CRF, 3CRF-D	Microconvex	JSM-115	Stainless Steel	Two Angles	11-23G
GE Healthcare	C358, 3C (3Cb), 3C-RS	Convex	JSM-157	Stainless Steel	Three Angles	11-23G
GE Healthcare	3S, 3S-RS, 3S-D, 5S, 5S-RS	Phased	JSM-163	Stainless Steel	Three Angles	11-23G
GE Healthcare	E72, E721, MTZ, MZ	Endocavity	JEM-041	Stainless Steel	Single Angle	16-18G
GE Healthcare	C4-8-D	Convex	JSM-172	Stainless Steel	Three Angles	11-23G
GE Healthcare	RAB4-8-D, RAB4-8-RS, RAB4-8L, RAB2-5-D, RAB2-5-RS, RAB2-5L, 4D3C-L	4D	JSM-174	Stainless Steel	Single Angle	11-23G
GE Healthcare	3Sp-D	Phased	JSM-179	Stainless Steel	Three Angles	11-23G
GE Healthcare	AB2-5, AB2-7, AB2-7-D, AB2-7-RS	Convex	JSM-180	Stainless Steel	Single Angle	11-23G
GE Healthcare	SP6-12	Linear	JSM-181	Stainless Steel	Single Angle	11-23G
GE Healthcare	C1-6, C1-6-D, C1-6VN-D	Convex	JSM-061	Stainless Steel	Three Angles	11-23G
GE Healthcare	C2-9-D, C2-9VN-D	Convex	JSM-182	Stainless Steel	Three Angles	11-23G

GE Healthcare	RAB6-D, RAB6-RS, RAB2-6-RS	4D	JSM-063	Stainless Steel	Three Angles	11-23G
GE Healthcare	E7C-RC	Endocavity	JEM-007	Stainless Steel	Single Angle	16-18G
GE Healthcare	C2-5	Convex	JSM-188	Stainless Steel	Three Angles	11-23G
GE Healthcare	IC9, IC9-RS	Endocavity	JEM-055	Stainless Steel	Single Angle	16-18G
GE Healthcare	ERB, E7C8L-RS	Transvaginal	JSM-129	Stainless Steel	Single Angle	18G
GE Healthcare	BE9C, BE9CS, BE9CS-RS, BE9CS-D	Rectal	JEM-014	Stainless Steel	Single Angle	16-18G
GE Healthcare	M5Sc, M5Sc-D, M5Sc-RS	Phased	JSM-136	Stainless Steel	Three Angles	11-23G
GE Healthcare	8C, 8C-RS	Microconvex	JSM-112	Stainless Steel	Free Angle	18-20G
GE Healthcare	M3S, M4S	Phased	JSM-131	Stainless Steel	Three Angles	11-23G
GE Healthcare	RM6C	4D	JSM-160	Stainless Steel	Single Angle	11-23G
GE Healthcare	3Sc, 3Sc-RS	Phased	JSM-130	Stainless Steel	Three Angles	11-23G
GE Healthcare	L8-18i-D, L8-18i-RS	Linear	JSM-111	Stainless Steel	Free Angle	18-20G
GE Healthcare	C2-7-D, C2-7VN-D	Microconvex	JSM-016	Stainless Steel	Two Angles	11-23G
GE Healthcare	C2-6b-D	Convex	JSM-104	Stainless Steel	Four Angles	18-19G
GE Healthcare	L3-12-D, L3-12-RS	Linear	JSM-039	Stainless Steel	Three Angles	11-23G
GE Healthcare	eM6C	4D	JSM-271	Stainless Steel	Three Angles	11-23G
GE Healthcare	L2-9-D, L2-9VN-D	Linear	JSM-327	Stainless Steel	Five Angles	11-23G
GE Healthcare	4Vc, 4Vc-D	Phased	JSM-328	Stainless Steel	Three Angles	11-23G
GE Healthcare	C3-10-D	Microconvex	JSM-122	Stainless Steel	Free Angle	11-23G
GE Healthcare	12S-RS, 12S-D	Phased	JSM-352	Stainless Steel	Free Angle	11-23G
GE Healthcare	L6-24-D	Linear	JSM-398	Stainless Steel	Free Angle	18-20G
GE Healthcare	6S-RS, 6S-D	Phased	JSM-433	Stainless Steel	Free Angle	11-23G
GE Healthcare	ML4-20, ML4-20-D, ML4-20VN-D	Linear	JSM-384	Stainless Steel	Five Angles	11-23G
Konica Minolta	EC9-3	Endocavity	JEM-010	Stainless Steel	Single Angle	16-18G
Konica Minolta	C5-2	Convex	JSM-056	Stainless Steel	Two Angles	11-23G
Konica Minolta	L18-4	Linear	JSM-057	Stainless Steel	Three Angles	14-22G
Konica Minolta	S4-2	Phased	JSM-021	Stainless Steel	Free Angle	11-23G
Konica Minolta	MC10-3	Microconvex	JSM-272	Stainless Steel	Two Angles	14-22G
Konica Minolta	L14-4	Linear	JSM-250	Stainless Steel	Two Angles	11-23G
Konica Minolta	L11-3	Linear	JSM-251	Stainless Steel	Three Angles	14-22G
Konica Minolta	HL18-4	Linear	JSM-252	Stainless Steel	Free Angle	18-20G

Konica Minolta	WL13-3	Linear	JSM-346	Stainless Steel	Three Angles	11-23G
Mindray-Zonare	E9-4	Endocavity	JEM-042	Stainless Steel	Single Angle	16-18G
Mindray-Zonare	C6-2	Convex	JSM-116	Stainless Steel	Single Angle	16G/18G
Mindray-Zonare	C9-3	Microconvex	JSM-229	Stainless Steel	Free Angle	16G/18G
Mindray	65EC10(EA,EB,EC,ED,HA), 65EB10EA, 6CV1(s,P), CB10- 4(s,P,E,A), V10-4(s), V10- 4B(s,P), V11-3(E,s), V11- 3B(E), V11-3W(s,E), E9-3, V11-3H(E,U,s)	Endocavity	JEM-028	Stainless Steel	Single Angle	16-18G
Mindray	C5-2, C5-2s, C5-2E	Convex	JSM-015	Stainless Steel	Three Angles	11-23G
Mindray	2P2(s), 2P2P, P4-2(s), P4-2E, P4-2NE, P4-2Ns, SP5-1s, SP5- 1E, SP5-1U, SP5-1N, SP5-1Ns	Phased	JSM-011	Stainless Steel	Two Angles	11-23G
Mindray	3C5A, 3C5, 3C5s, 3C5P	Convex	JSM-006	Stainless Steel	Three Angles	11-23G
Mindray	7L4(s,A,P), 7L4B(P), 10L4(s), 7L5(s,P), L7-3(s,E), L11-4(s), L12-4(s), L12-3E, L13-3(s), L14-6N(s,E), L14-6W(s,E,U), 75L53EA(HA)	Linear	JSM-007	Stainless Steel	Three Angles	11-23G
Mindray	35C50EA, 35C50EB, 35C50P	Convex	JSM-001	Stainless Steel	Three Angles	11-23G
Mindray	3C1, 3C1s	Microconvex	JSM-074	Stainless Steel	Two Angles	11-23G
Mindray	65C15EA, 6C2, 6C2s, 6C2P	Microconvex	JSM-075	Stainless Steel	Two Angles	16G/18G
Mindray	DE10-3(s,E,U), DE10- 3W(s,E,U)	Endocavity	JEM-063	Stainless Steel	Single Angle	16-18G
Mindray	75L38EA, 75L38EB, 75L38P	Linear	JSM-088	Stainless Steel	Three Angles	11-23G
Mindray	LM15-6U	Linear	JSM-036	Stainless Steel	Three Angles	11-23G
Mindray	L14-5W, L14-5Ws, L14-5WE, L14-5WU	Linear	JSM-127	Stainless Steel	Four Angles	11-23G
Mindray	C4-1, C4-1U, C4-1s	Microconvex	JSM-123	Stainless Steel	Three Angles	14-22G
Mindray	LM14-6E, LM14-6s, LM16-4U	Linear	JSM-237	Stainless Steel	Three Angles	11-23G
Mindray	SC5-1E	Convex	JSM-045	Stainless Steel	Three Angles	11-23G
Mindray	C5-1(s,E,U), SC6-1(s,E,U), C6- 2(s,P), SC5-1N(s), SC5-1s	Convex	JSM-022	Stainless Steel	Three Angles	11-23G
Mindray	L10-3E, L10-3s	Linear	JSM-284	Stainless Steel	Three Angles	11-23G
Mindray	C7-3E	Convex	JSM-219	Stainless Steel	Three Angles	11-23G
Mindray	C11-3(s, E, U)	Microconvex	JSM-118	Stainless Steel	Three Angles	11-23G
Mindray	SC5-1U, SC5-1E	Convex	JSM-031	Stainless Steel	Three Angles	11-23G
Mindray	D8-2E, D8-2U	4D	JSM-281	Stainless Steel	Three Angles	14-22G
Mindray	D8-4U	4D	JSM-144	Stainless Steel	Single Angle	16G/18G
Mindray	C6-2GE, C6-2Gs, C6-2GU	Microconvex	JSM-024	Stainless Steel	Three Angles	11-23G

Mindray	P7-3, P7-3s, P7-3E, P7-3P, P7-3U	Phased	JSM-221	Stainless Steel	Two Angles	11-23G
Mindray	L14-3U, L13-3Ns	Linear	JSM-214	Stainless Steel	Three Angles	11-23G
Mindray	L11-3U	Linear	JSM-126	Stainless Steel	Three Angles	11-23G
Mindray	L20-5, L20-5s, L20-5E, L20-5U	Linear	JSM-287	Stainless Steel	Three Angles	11-23G
Mindray	V11-3HB	Endocavity	JEM-043	Stainless Steel	Single Angle	16-18G
Mindray	6LB7(s), 6LB7E, 6LB7P, 6LE7(s), 6LE7P, 65EL60EA	Transvaginal	JSM-109	Stainless Steel	Single Angle	18G
Mindray	L16-4Hs, L16-4HE, L16-4HU	Linear	JSM-204	Stainless Steel	Single Angle	16G/18G
Mindray	L9-3U, L9-3E, L9-3s	Linear	JSM-264	Stainless Steel	Three Angles	11-23G
Mindray	SC8-2U, SC8-2s	Convex	JSM-265	Stainless Steel	Three Angles	11-23G
Mindray	7LT4(s), 7LT4P, 75LT38EA, 75LT40HA	Linear	JSM-282	Stainless Steel	Single Angle	18G/20G
Mindray	L14-6(s), L14-6P, L14-6E, 10L24EA	Linear	JSM-274	Stainless Steel	Three Angles	14-22G
Mindray	L11-3VNs	Linear	JSM-074	Stainless Steel	Three Angles	11-23G
Mindray	L12-3RC, L12-3RCs	Linear	JSM-146	Stainless Steel	Free Angle	20-23G
Mindray	ELC13-4U, ELC13-4s	Transvaginal	JSM-246	Stainless Steel	Single Angle	18G
Mindray	L14-3WU, L14-3WE, L14-3Ws, L15-3WU, L13-3WE, L13-	Linear	JSM-305	Stainless Steel	Three Angles	11-23G
Mindray	V11-3H(E,U,s)	Endocavity	JEM-064	Stainless Steel	Single Angle	16-18G
Mindray	DE11-3U, DE11-3P, DE11-3s, DE11-3E, DE11-3	Endocavity	JEM-008	Stainless Steel	Single Angle	16-18G
Mindray	75L60HB, 75L60EA, 7L6(s)	Linear	JSM-415	Stainless Steel	Three Angles	11-23G
Mindray	35C20EA	Microconvex	JSM-298	Stainless Steel	Two Angles	11-23G
Mindray	C6-1	Convex	JSM-348	Stainless Steel	Three Angles	11-23G
Mindray	SCM7-1U	Convex	JSM-444	Stainless Steel	Three Angles	11-23G
Mindray	SL10-3U	Linear	JSM-446	Stainless Steel	Three Angles	11-23G
Mindray	LM18-4WU	Linear	JSM-447	Stainless Steel	Three Angles	11-23G
Mindray	SPM6-1U	Phased	JSM-443	Stainless Steel	Two Angles	11-23G
Mindray	SC9-2U	Convex	JSM-445	Stainless Steel	Three Angles	11-23G
Mindray	SC6-1GU	Microconvex	JSM-448	Stainless Steel	Three Angles	11-23G
Mindray	SV10-2U	Endocavity	JEM-092	Stainless Steel	Single Angle	16-18G
Mindray	DE11-3WU, DE11-3WE, DE11-3Ws	Endocavity	JEM-080	Stainless Steel	Single Angle	16-18G
Mindray	SDM9-2U	4D	JSM-451	Stainless Steel	Three Angles	11-23G
Mindray	SDE10-2WU	Endocavity	JEM-083	Stainless Steel	Single Angle	16-18G
Neusoft	SC1-6, SC1-6H	Convex	JSM-159	Stainless Steel	Three Angles	11-23G

Neusoft	L8-17X	Linear	JSM-199	Stainless Steel	Three Angles	11-23G
Neusoft	E3-10	Endocavity	JEM-030	Stainless Steel	Single Angle	16-18G
Neusoft	EV3-10X, EC3-10X	Endocavity	JEM-033	Stainless Steel	Single Angle	16-18G
Neusoft	C5-2	Convex	JSM-288	Stainless Steel	Three Angles	14-23G
Neusoft	EC8-4	Endocavity	JEM-015	Stainless Steel	Single Angle	16-18G
Neusoft	L12-5, L17-7H	Linear	JSM-240	Stainless Steel	Three Angles	14-23G
Neusoft	EV10-3H	Endocavity	JEM-066	Stainless Steel	Single Angle	16-18G
Neusoft	EC10-3H	Endocavity	JEM-067	Stainless Steel	Single Angle	16-18G
Neusoft	C10-3cc	Endocavity	JEM-065	Stainless Steel	Single Angle	16-18G
Neusoft	C10-3cc	Endocavity	JEM-065	Stainless Steel	Single Angle	16-18G
Neusoft	C8-4	Microconvex	JSM-263	Stainless Steel	Three Angles	14-23G
Neusoft	S4-2M	Phased	JSM-337	Stainless Steel	Two Angles	11-23G
Neusoft	L12-3M	Linear	JSM-338	Stainless Steel	Three Angles	11-23G
Neusoft	ECL12-3	Transvaginal	JSM-347	Stainless Steel	Single Angle	18G
Neusoft	SC6-1H	Convex	JSM-414	Stainless Steel	Four Angles	11-23G
Philips	C8-4v, C9-3v, C10-3v	Endocavity	JEM-023	Stainless Steel	Single Angle	16-18G
Philips	L12-5 50	Linear	JSM-046	Stainless Steel	Five Angles	11-23G
Philips	C5-2(HD5, HD7, HD11, HD11XE, HDI 4000, HDI 5000/5500cv, iE33 and iU22)	Convex	JSM-073	Stainless Steel	Single Angle	11-23G
Philips	C5-2(ENVISOR, HD6, HD7XE, HD9 and HD15)	Convex	JSM-073	Stainless Steel	Four Angles	11-23G
Philips	C5-1(CX50, EPIQ, Affiniti)	Convex	JSM-069	Stainless Steel	Four Angles	11-23G
Philips	C5-1(IU22, IE33, HD15)	Convex	JSM-069	Stainless Steel	Four Angles	11-23G
Philips	C5-1(For Verza)	Convex	JSM-069	Stainless Steel	Five Angles	11-23G
Philips	C5-2(ClearVue)	Convex	JSM-087	Stainless Steel	Four Angles	11-23G
Philips	PA4-2, S4-2	Phased	JSM-026	Stainless Steel	Two Angles	11-23G
Philips	C9-4v(Affiniti, SPARQ, InnoSight, 3300)	Endocavity	JEM-021	Stainless Steel	Single Angle	16-18G
Philips	S4-1, P4-1	Phased	JSM-094	Stainless Steel	Two Angles	11-23G
Philips	C9-4v(ClearVue)	Endocavity	JEM-035	Stainless Steel	Single Angle	16-18G
Philips	L12-4(Affiniti and ClearVue), L12-3(Affiniti, EPIQ and CX50)	Linear	JSM-044	Stainless Steel	Three Angles	11-23G
Philips	S4-1(ClearVue)	Phased	JSM-096	Stainless Steel	Three Angles	11-23G
Philips	L12-3(HD5, HD6, HD7, HD11, HD15 and HD7XE)	Linear	JSM-099	Stainless Steel	Single Angle	11-23G
Philips	S5-1	Phased	JSM-034	Stainless Steel	Three Angles	11-23G

Philips	L9-3	Linear	JSM-153	Stainless Steel	Two Angles	11-23G
Philips	C9-5ec, C9-5 ICT	Rectal	JEM-001	Stainless Steel	Single Angle	16-18G
Philips	L11-3	Linear	JSM-148	Stainless Steel	Free Angle	18-20G
Philips	C3540	Convex	JSM-168	Stainless Steel	Single Angle	11-23G
Philips	L17-5, L18-5	Linear	JSM-147	Stainless Steel	Three Angles	11-23G
Philips	L12-4(CX30, SPARQ, InnoSight)	Linear	JSM-170	Stainless Steel	Three Angles	11-23G
Philips	C6-2, C6-3	Convex	JSM-173	Stainless Steel	Four Angles	11-23G
Philips	C8-5(CX30, CX50, HD9 and HD15)	Microconvex	JSM-175	Stainless Steel	Three Angles	11-23G
Philips	C8-5(Affiniti 50, Affiniti 70, EPIQ 5, EPIQ 7, ENVISOR, HD7XE, HD11, IU22, IE33)	Microconvex	JSM-175	Stainless Steel	Single Angle	11-23G
Philips	X6-1	Matrix	JSM-154	Stainless Steel	Three Angles	11-23G
Philips	3D9-3v	Endocavity	JEM-040	Stainless Steel	Single Angle	16-18G
Philips	VL13-5	4D	JSM-238	Stainless Steel	Two Angles	11-23G
Philips	C9-4ec(HD3)	Endocavity	JEM-022	Stainless Steel	Single Angle	16-18G
Philips	L12-4(HD8)	Linear	JSM-318	Stainless Steel	Single Angle	11-23G
Philips	C10-4ec	Rectal	JEM-004	Stainless Steel	Single Angle	16-18G
Philips	C9-2	Convex	JSM-162	Stainless Steel	Four Angles	11-23G
Philips	BP10-5ec, C9-4ec(HD7 XE & HD9)	Endocavity	JEM-030	Stainless Steel	Single Angle	16-18G
Philips	X5-1	Matrix	JSM-292	Stainless Steel	Free Angle	14-22G
Philips	V6-2	4D	JSM-210	Stainless Steel	Two Angles	11-23G
Philips	S8-3	Phased	JSM-120	Stainless Steel	Free Angle	11-23G
Philips	C5-2(VISIQ)	Convex	JSM-257	Stainless Steel	Free Angle	11-23G
Philips	eL18-4	Linear	JSM-084	Stainless Steel	Five Angles	11-23G
Philips	XL14-3	Matrix	JSM-358	Stainless Steel	Five Angles	11-23G
Philips	V9-2	4D	JSM-360	Stainless Steel	Five Angles	14-22G
Philips	L12-5 38(Affiniti 50, ClearVue)	Linear	JSM-145	Stainless Steel	Free Angle	18-20G
Philips	C5-2(Lumify)	Convex	JSM-311	Stainless Steel	Free Angle	11-23G
Philips	L12-4(Lumify)	Linear	JSM-312	Stainless Steel	Free Angle	16-18G
Philips	S4-1(Lumify)	Phased	JSM-313	Stainless Steel	Free Angle	14-22G
Philips	mC7-2	Microconvex	JSM-335	Stainless Steel	Five Angles	16G/18G
Philips	L15-7io	Linear	JSM-252	Stainless Steel	Free Angle	22-23G
Philips	L12-3 Ergo	Linear	JSM-366	Stainless Steel	Free Angle	11-23G

Philips	S9-2	Phased	JSM-494	Stainless Steel	Free Angle	11-23G
Samsung	VR5-9	Endocavity	JEM-037	Stainless Steel	Single Angle	16-18G
Samsung	CS1-4	Convex	JSM-003	Stainless Steel	Three Angles	11-23G
Samsung	C2-6, C2-6IC	Convex	JSM-062	Stainless Steel	Three Angles	11-23G
Samsung	V5-9, EV3-10B	Endocavity	JEM-063	Stainless Steel	Single Angle	16-18G
Samsung	C1-4EC	Convex	JSM-019	Stainless Steel	Two Angles	11-23G
Samsung	L3-8, LN5-12/40	Linear	JSM-004	Stainless Steel	Single Angle	14-22G
Samsung	C2-8, CN2-8, C2-8-N, C2-8N	Convex	JSM-083	Stainless Steel	Single Angle	11-23G
Samsung	P1-4, PE2-4	Phased	JSM-026	Stainless Steel	Two Angles	11-23G
Samsung	L5-13, L5-13IS, L7-16IS, L7-16	Linear	JSM-103	Stainless Steel	Three Angles	11-23G
Samsung	ER4-9/10ED, EV4-9/10ED, EVN4-9(-H), EV4-9, ER4-9,	Endocavity	JEM-027	Stainless Steel	Single Angle	16-18G
Samsung	LM4-15B, LM4-15BD	Linear	JSM-036	Stainless Steel	Three Angles	11-23G
Samsung	HL5-12ED	Linear	JSM-110	Stainless Steel	Three Angles	11-23G
Samsung	EC4-9IS, EC4-9ES, ER4-9ES	Endocavity	JEM-011	Stainless Steel	Single Angle	16-18G
Samsung	EA2-11B	Endocavity	JEM-060	Stainless Steel	Single Angle	16-18G
Samsung	9L, 9L-D, 9L-RS	Linear	JSM-048	Stainless Steel	Three Angles	11-23G
Samsung	3D4-9, 3D4-9ES	Endocavity	JEM-054	Stainless Steel	Single Angle	16-18G
Samsung	EC4-9/10ED, EC4-9/13CD, EC4-9/10R, EC4-9ED	Endocavity	JEM-022	Stainless Steel	Single Angle	16-18G
Samsung	CF4-9	Microconvex	JSM-029	Stainless Steel	Two Angles	11-23G
Samsung	LN5-12, LE5-12	Linear	JSM-149	Stainless Steel	Three Angles	11-23G
Samsung	C3-7ED, C3-7IM, C3-7EP	Convex	JSM-289	Stainless Steel	Three Angles	11-23G
Samsung	V4-8, VE4-8	4D	JSM-144	Stainless Steel	Single Angle	16G/18G
Samsung	NEV4-9ES, NER4-9ES, E3-12A	Endocavity	JEM-030	Stainless Steel	Single Angle	16-18G
Samsung	CA1-7A, CA1-7AD	Convex	JSM-071	Stainless Steel	Four Angles	11-23G
Samsung	L3-12A, LA3-12A, LA3-14AD, L5-13/50, LA2-14A	Linear	JSM-072	Stainless Steel	Three Angles	11-23G
Samsung	P3-8CA, PA3-8B	Phased	JSM-120	Stainless Steel	Free Angle	11-23G
Samsung	PA1-5A	Phased	JSM-092	Stainless Steel	Free Angle	14-22G
Samsung	SC1-6	Convex	JSM-217	Stainless Steel	Two Angles	11-23G
Samsung	C2-8A, CA2-8A, CA2-8AD, CA2-9A, CA2-9AD, CF2-8	Convex	JSM-091	Stainless Steel	Three Angles	11-23G
Samsung	CA3-10A	Convex	JSM-137	Stainless Steel	Four Angles	11-23G
Samsung	LA4-18B, LA4-18BD	Linear	JSM-245	Stainless Steel	Three Angles	11-23G

Samsung	CV1-8A, CV1-8AD	4D	JSM-161	Stainless Steel	Three Angles	14-22G
Samsung	LA3-16A, LA3-16AD, LF5-13	Linear	JSM-090	Stainless Steel	Three Angles	11-23G
Samsung	C2-5	Convex	JSM-315	Stainless Steel	Three Angles	11-23G
Samsung	EA2-11AR, EA2-11AV, EA2-11ARD, EA2-11AVD, EA2-11ARE, EA2-11AVE	Endocavity	JEM-070	Stainless Steel	Single Angle	16-18G
Samsung	EA2-11AR, EA2-11AV, EA2-11ARD, EA2-11AVD, EA2-11ARE, EA2-11AVE	Endocavity	JEM-070	Stainless Steel	Single Angle	16-18G
Samsung	EV2-10A	Endocavity	JEM-073	Stainless Steel	Single Angle	16-18G
Samsung	EV2-10A	Endocavity	JEM-073	Stainless Steel	Single Angle	16-18G
Samsung	CA1-7S, CA1-7SD	Convex	JSM-389	Stainless Steel	Five Angles	11-23G
Samsung	miniER7	Endocavity	JEM-081	Stainless Steel	Single Angle	16-18G
Samsung	LA4-18A, LA4-18AD	Linear	JSM-322	Stainless Steel	Five Angles	11-23G
Siemens	EC9-4(ACUSON S1000, S2000, S3000, Antares)	Endocavity	JEM-018	Stainless Steel	Single Angle	16-18G
Siemens	CH6-2, 5C2g, 6C2, 7C2	Convex	JSM-062	Stainless Steel	Single Angle	11-23G
Siemens	VF13-5	Linear	JSM-078	Stainless Steel	Single Angle	11-23G
Siemens	EC9-4, EV9-4(ACUSON X150, X300, P500, OMNI, NX2, NX3; SONOLINE G40)	Endocavity	JEM-029	Stainless Steel	Single Angle	16-18G
Siemens	4C1, 6C1 HD	Convex	JSM-019	Stainless Steel	Two Angles	11-23G
Siemens	CH4-1, CH5-2, 5C1(Juniper)	Convex	JSM-081	Stainless Steel	Two Angles	11-23G
Siemens	9L4, VFX9-4, VF7-3, VF10-5, 11L4	Linear	JSM-008	Stainless Steel	Two Angles	11-23G
Siemens	4V1, PH4-1	Phased	JSM-043	Stainless Steel	Two Angles	11-23G
Siemens	6L3, 8L5	Linear	JSM-060	Stainless Steel	Single Angle	11-23G
Siemens	EC-10C5, EC9-4, EV9-4(ACUSON X500; SONOLINE	Endocavity	JEM-011	Stainless Steel	Single Angle	16-18G
Siemens	9VE3, EV9F3	Endocavity	JEM-054	Stainless Steel	Single Angle	16-18G
Siemens	4V1c	Phased	JSM-164	Stainless Steel	Free Angle	14-22G
Siemens	10L4	Linear	JSM-294	Stainless Steel	Five Angles	11-23G
Siemens	5C1(Sequoia), 6C1(Juniper)	Convex	JSM-291	Stainless Steel	Five Angles	11-23G
Siemens	14L5, VFX13-5	Linear	JSM-191	Stainless Steel	Single Angle	11-23G
Siemens	EV-8C4, EVC8	T	JEM-039	Stainless Steel	Single Angle	16-18G
Siemens	DAX	Matrix	JSM-293	Stainless Steel	Five Angles	11-23G
Siemens	MC9-4, EC9-4W, 9MC3, 9EC4	Endocavity	JEM-010	Stainless Steel	Single Angle	16-18G
Siemens	18L6HD, 18L6	Linear	JSM-058	Stainless Steel	Two Angles	11-23G
Siemens	18L6HD, 18L6(OXANA)	Linear	JSM-058	Stainless Steel	Two Angles	11-23G

Siemens	10MC3	Endocavity	JEM-052	Stainless Steel	Single Angle	16-18G
Siemens	9EVF4, EV9F4	Endocavity	JEM-045	Stainless Steel	Single Angle	16-18G
Siemens	8C3HD, 9C3	Convex	JSM-239	Stainless Steel	Two Angles	11-23G
Siemens	VF12-4, 12L4, 12L3	Linear	JSM-142	Stainless Steel	Single Angle	11-23G
Siemens	BP10-3	Endocavity	JEM-049	Stainless Steel	Single Angle	16-18G
Siemens	L13-5(freestyle)	Linear	JSM-244	Stainless Steel	Free Angle	11-23G
Siemens	EV10-5, EC10-5w	Endocavity	JEM-026	Stainless Steel	Single Angle	16-18G
Siemens	BP9-4	Rectal	JEM-062	Stainless Steel	Single Angle	16-18G
Siemens	8V3	Phased	JSM-119	Stainless Steel	Free Angle	14-22G
Siemens	11M3	Microconvex	JSM-272	Stainless Steel	Two Angles	14-22G
Siemens	14L5(Sequoia)	Linear	JSM-373	Stainless Steel	Five Angles	11-23G
Siemens	5V1	Phased	JSM-363	Stainless Steel	Free Angle	11-23G
Siemens	VF16-5, 16L4	Linear	JSM-468	Stainless Steel	Free Angle	18-20G
SIUI	V6L11C, V6L11K, V6I11F, V5H11, V6LA, V6LC	Endocavity	JEM-002	Stainless Steel	Single Angle	16-18G
SIUI	C3L60, C3LC, C3L60C	Convex	JSM-150	Stainless Steel	Three Angles	11-23G
SIUI	L8L50	Linear	JSM-141	Stainless Steel	Three Angles	14-22G
SIUI	L8L38, L10LC	Linear	JSM-319	Stainless Steel	Three Angles	14-22G
SIUI	L8LC(50), L8LN(50)	Linear	JSM-222	Stainless Steel	Three Angles	11-23G
SIUI	C3LC, C3L60C	Convex	JSM-336	Stainless Steel	Three Angles	11-23G
SIUI	L8L38C	Linear	JSM-278	Stainless Steel	Free Angle	18-20G
SIUI	V6LC (4D)	Endocavity	JEM-090	Stainless Steel	Single Angle	16-18G
SIUI	C3L60	Convex	JSM-336	Stainless Steel	Three Angles	11-23G
SonoScape	6V1	Endocavity	JEM-017	Stainless Steel	Single Angle	16-18G
SonoScape	6V3	Endocavity	JEM-003	Stainless Steel	Single Angle	16-18G
SonoScape	C322	Microconvex	JSM-166	Stainless Steel	Two Angles	11-23G
SonoScape	C351	Convex	JSM-183	Stainless Steel	Three Angles	11-23G
SonoScape	6V1	Endocavity	JEM-006	Stainless Steel	Single Angle	16-18G
SonoScape	6V3	Endocavity	JEM-031	Stainless Steel	Single Angle	16-18G
SonoScape	6V4	Endocavity	JEM-009	Stainless Steel	Single Angle	16-18G
SonoScape	C344	Convex	JSM-185	Stainless Steel	Single Angle	11-23G
SonoScape	C612	Microconvex	JSM-186	Stainless Steel	Two Angles	11-23G

SonoScape	6V7	Endocavity	JEM-047	Stainless Steel	Single Angle	16-18G
SonoScape	C613	Microconvex	JSM-189	Stainless Steel	Three Angles	11-23G
SonoScape	3C-A	Convex	JSM-190	Stainless Steel	Four Angles	11-23G
SonoScape	4P-A, 3P-A	Phased	JSM-135	Stainless Steel	Two Angles	11-23G
SonoScape	BCL10-5	Transvaginal	JSM-330	Stainless Steel	Single Angle	18G
SonoScape	L742(II), 10L1(II)	Linear	JSM-207	Stainless Steel	Two Angles	11-23G
SonoScape	L741, L743, L746	Linear	JSM-203	Stainless Steel	Single Angle	11-23G
SonoScape	7P-A, 7P-B, 5P2	Phased	JSM-205	Stainless Steel	Three Angles	14-22G
SonoScape	S1-5	Phased	JSM-231	Stainless Steel	Three Angles	14-22G
SonoScape	C1-6	Convex	JSM-206	Stainless Steel	Four Angles	11-23G
SonoScape	L742	Linear	JSM-270	Stainless Steel	Two Angles	11-23G
SonoScape	C353	Convex	JSM-208	Stainless Steel	Two Angles	11-23G
SonoScape	VC2-9	4D	JSM-281	Stainless Steel	Three Angles	14-22G
SonoScape	10L-I	Linear	JSM-277	Stainless Steel	Three Angles	11-23G
SonoScape	12L-A	Linear	JSM-216	Stainless Steel	Three Angles	11-23G
SonoScape	VE9-5, VE3-10	Endocavity	JEM-046	Stainless Steel	Single Angle	16-18G
SonoScape	MC1-6	Microconvex	JSM-407	Stainless Steel	Two Angles	11-23G
SonoScape	BCC9-5, BCC9-4	Endocavity	JEM-053	Stainless Steel	Single Angle	16-18G
SonoScape	12L-B	Linear	JSM-275	Stainless Steel	Three Angles	11-23G
SonoScape	C1-6A	Convex	JSM-306	Stainless Steel	Four Angles	11-23G
SonoScape	L752	Linear	JSM-152	Stainless Steel	Single Angle	11-23G
SonoScape	9L-A	Linear	JSM-171	Stainless Steel	Three Angles	11-23G
SonoScape	EC9-5	Rectal	JEM-072	Stainless Steel	Single Angle	16-18G
SonoScape	13L-A	Linear	JSM-331	Stainless Steel	Three Angles	11-23G
SonoScape	18L-A	Linear	JSM-332	Stainless Steel	Three Angles	11-23G
SonoScape	C2-9	Convex	JSM-333	Stainless Steel	Four Angles	11-23G
SonoScape	VC6-2	Convex	JSM-372	Stainless Steel	Three Angles	11-23G
SonoScape	6V5, EC2	Endocavity	JEM-087	Stainless Steel	Single Angle	16-18G
SonoScape	C361, C362	Convex	JSM-357	Stainless Steel	Single Angle	11-23G
SonoScape	C1-5	Convex	JSM-408	Stainless Steel	Three Angles	11-23G
SonoScape	L2-9	Linear	JSM-409	Stainless Steel	Three Angles	11-23G

SonoScape	L3-11	Linear	JSM-410	Stainless Steel	Three Angles	11-23G
SonoScape	VC2-9A	4D	JSM-472	Stainless Steel	Three Angles	11-23G
SuperSonic	SC6-1	Convex	JSM-176	Stainless Steel	Four Angles	11-23G
SuperSonic	SL15-4	Linear	JSM-066	Stainless Steel	Free Angle	18-20G
SuperSonic	SE12-3	Endocavity	JEM-025	Stainless Steel	Single Angle	16-18G
SuperSonic	SL10-2(I)	Linear	JSM-038	Stainless Steel	Free Angle	18-20G
SuperSonic	SL10-2(II)	Linear	JSM-038	Stainless Steel	Free Angle	11-23G
SuperSonic	XC6-1, SC6-1(II), C6-1X	Convex	JSM-158	Stainless Steel	Four Angles	11-23G
SuperSonic	SMC12-3	Microconvex	JSM-023	Stainless Steel	Free Angle	14-22G
SuperSonic	C9-2X	Convex	JSM-262	Stainless Steel	Four Angles	11-23G
SuperSonic	XP5-1	Phased	JSM-139	Stainless Steel	Free Angle	14-22G
SuperSonic	SL15-4(II), SL18-5	Linear	JSM-268	Stainless Steel	Free Angle	18-20G
SuperSonic	SMC12-3(II)	Microconvex	JSM-023	Stainless Steel	Free Angle	14-22G
Terason	5C2A	Convex	JSM-089	Stainless Steel	Single Angle	11-23G
Terason	4V2, 4V2A	Phased	JSM-026	Stainless Steel	Two Angles	11-23G
Terason	8EC4A	Endocavity	JEM-036	Stainless Steel	Single Angle	16-18G
Terason	12L5-A	Linear	JSM-105	Stainless Steel	Single Angle	11-23G
Terason	16L5	Linear	JSM-297	Stainless Steel	Three Angles	11-23G
Terason	15L4A	Linear	JSM-318	Stainless Steel	Single Angle	11-23G
Terason	8BP4	Transvaginal	JSM-013	Stainless Steel	Single Angle	14/16G, 18/20G
Terason	8V3, 8V3A	Phased	JSM-120	Stainless Steel	Free Angle	11-23G
VINNO	G4-9E	Endocavity	JEM-028	Stainless Steel	Single Angle	16-18G
VINNO	F4-12L, A4-12L	Linear	JSM-093	Stainless Steel	Three Angles	11-23G
VINNO	F2-5C, F2-5CE, A2-5C	Convex	JSM-050	Stainless Steel	Three Angles	11-23G
VINNO	BP4-9	Transvaginal	JSM-013	Stainless Steel	Single Angle	18G
VINNO	X6-16L	Linear	JSM-464	Stainless Steel	Three Angles	11-23G
VINNO	S1-8C	Convex	JSM-259	Stainless Steel	Three Angles	11-23G
VINNO	G2-5C	Convex	JSM-223	Stainless Steel	Three Angles	11-23G
VINNO	X4-12L	Linear	JSM-260	Stainless Steel	Three Angles	11-23G
VINNO	X4-9E	Endocavity	JEM-048	Stainless Steel	Single Angle	16-18G
VINNO	F4-9E, A4-9E	Endocavity	JEM-099	Stainless Steel	Single Angle	16-18G

VINNO	S2-9C	Convex	JSM-499	Stainless Steel	Three Angles	11-23G
VINNO	U5-15L, U5-15LE	Linear	JSM-435	Stainless Steel	Free Angle	11-23G
VINNO	G3-8M	Microconvex	JSM-483	Stainless Steel	Two Angles	11-23G
VINNO	I7-18L	Linear	JSM-453	Stainless Steel	Free Angle	16-18G
VINNO	S1-6P(V)	Phased	JSM-485	Stainless Steel	Three Angles	11-23G
Wisonic	S4-2	Phased	JSM-026	Stainless Steel	Two Angles	11-23G
Wisonic	C6-1B-H	Convex	JSM-295	Stainless Steel	Three Angles	11-23G
Wisonic	L15-4WB-H	Linear	JSM-296	Stainless Steel	Three Angles	11-23G
Wisonic	EV10-4, EV11-3-H	Endocavity	JEM-028	Stainless Steel	Single Angle	16-18G
Wisonic	L15-6-H	Linear	JSM-297	Stainless Steel	Three Angles	11-23G
Wisonic	C7-3B, C5-1B, C5-2B	Convex	JSM-059	Stainless Steel	Three Angles	11-23G
Wisonic	L15-4B, L15-4NB	Linear	JSM-040	Stainless Steel	Three Angles	11-23G
Wisonic	L15-4	Linear	JSM-041	Stainless Steel	Three Angles	11-23G
Wisonic	C5-1, C5-2, C7-3	Convex	JSM-050	Stainless Steel	Three Angles	11-23G
Wisonic	L10-5	Linear	JSM-032	Stainless Steel	Three Angles	11-23G
Wisonic	P4-1	Phased	JSM-020	Stainless Steel	Two Angles	11-23G
Wisonic	EV11-2-U	Endocavity	JEM-071	Stainless Steel	Single Angle	16-18G
Wisonic	C11-3-U	Microconvex	JSM-420	Stainless Steel	Two Angles	11-23G
Wisonic	L15-4W-U	Linear	JSM-426	Stainless Steel	Three Angles	11-23G
Wisonic	CS6-2-U	Convex	JSM-427	Stainless Steel	Three Angles	11-23G
Wisonic	C9-2G-U	Convex	JSM-428	Stainless Steel	Three Angles	11-23G
Wisonic	HU-11C3	Microconvex	JSM-452	Stainless Steel	Three Angles	11-23G
Wisonic	C8-3	Microconvex	JSM-483	Stainless Steel	Two Angles	11-23G
Wisonic	LH15-6-H	Linear	JSM-453	Stainless Steel	Free Angle	16-18G